

pH600





pH Economical Pocket Tester

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

pH600 pH Economical Pocket Tester pH600-AQ pH Economical Pocket Tester in Blister packaging

Milwaukee's economical testers are easy-to-use and low cost instruments to measure quick and reliable pH, EC or TDS values. Measuring electrical conductivity is the best way of checking the amount of salt or dissolved solids (TDS) in water.

Milwaukee provides you with a range of pocket testers that will allow you to measure from very low to very high conductivity solutions. All EC/TDS testers compensate for the temperature variance automatically.

The **pH600AQ** is a pH tester specifically designed for the aquarium market. The calibration procedure is very simple and fast. This tester has 1 point manual calibration, ± 0.1 pH accuracy and a range of 0.0 to 14.0 pH. The pH600 is supplied in a blister packaging complete with calibration screwdriver, batteries and instructions.

Specification	pH600
Range	0.0 to 14.0 pH
Resolution	0.1 pH
Accuracy	± 0.1 pH
Calibration	manual, 1 point
Environment	0 to 50°C / 32 to 122°F; max RH 95%
Battery Type / Battery Life	3 x 1.5V, alkaline / 700 hours of use
Packaging dimensions / weight	150 x 65 x 32 mm / 120 g

Accessories:

- **M10004B** pH 4.01 buffer solution 20 ml sachet (25 pcs)
- **M10007B** pH 7.01 buffer solution 20 ml sachet (25 pcs)
- **M10010B** pH 10.01 buffer solution 20 ml sachet (25 pcs)
- **MA9015** Electrode storage solution, 230 ml
- **MA9016** Electrode cleaning solution, 230 ml
- **M10000B** Electrode rinse solution, 20 mL sachet (25 pcs)
- **MA9305** 1.5V batteries (100 pcs)

Ordering information:

The pH600 is supplied in a blister packaging complete with calibration screwdriver, batteries and instructions.